


Product Datasheet

Actin antibody (orb1151844)

Description	Mouse monoclonal antibody to Actin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein was used as the immunogen for the Muscle Actin antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the Muscle Actin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the Muscle Actin antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Purity	Protein A/G affinity
Uniprot ID	P62736 , P68032 , P68133
Hazard Information	This Muscle Actin antibody is available for research use only.
Dilution Range	Western blot: 0.5-1ug/ml, Immunohistochemistry (FFPE): 0.5-2ug/ml
Expiration Date	12 months from date of receipt.



Western blot testing of human skeletal m...



SDS-PAGE analysis of purified, BSA-free ...