


Product Datasheet

HPV16 E2 antibody (orb1151847)

Description	Mouse monoclonal antibody to HPV16 E2
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Raised against Vaccinia-E2 followed by intravenous injection of the maltose binding protein MBP-E2
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the HPV16 E2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the HPV16 E2 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein A/G affinity
Uniprot ID	Not Applicable
Hazard Information	This HPV16 E2 antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml
Expiration Date	12 months from date of receipt.



SDS-PAGE analysis of purified, BSA-free ...