

Product Datasheet

Goat Anti-Rabbit IgG H&L/HRP antibody (orb1151906)

| | |
|----------------------------|---|
| Description | Goat Anti-Rabbit IgG H&L/HRP antibody |
| Species/Host | Goat |
| Reactivity | Rabbit |
| Conjugation | HRP |
| Tested Applications | WB |
| Immunogen | Native Rabbit IgG |
| Preservatives | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1.0 mg/ml |
| Storage | Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Dilution Range | WB=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |