

Product Datasheet

Histone H2B (Acetyl K11) antibody (orb1152019)

Description	Histone H2B (Acetyl K11) antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Histone H2B (Acetyl K11)
Target	Histone H2B (Acetyl K11)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	14 kDa
Uniprot ID	P33778
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.