



Product Datasheet

Catalase antibody (orb1152050)


Description	Catalase antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Catalase
Target	Catalase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1/Kappa
Clonality	Recombinant
MW	60 kDa
Uniprot ID	P04040
Dilution Range	WB=1:500-2000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P analysis of Human liver using Catalase antibody.



IHC-P analysis of rat kidney using Catalase antibody.



IHC-P analysis of mouse kidney using Catalase antibody.