

## Product Datasheet

### phospho-Tau (Thr217) antibody (orb1152153)

<b>Description</b>	phospho-Tau (Thr217) antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr217 (LP(p-T)PP)
<b>Target</b>	phospho-Tau (Thr217)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	79 kDa
<b>Uniprot ID</b>	<b>P10636</b>
<b>Expiration Date</b>	12 months from date of receipt.