

---

## Product Datasheet

**SARS-CoV-2 (2019-nCoV) Nucleocapsid protein,  
Biotin conjugated antibody (orb1152154)**

**Description**

SARS-CoV-2 (2019-nCoV)

	conjugated antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant SARS-CoV-2 Nucleocapsid Protei
<b>Target</b>	SARS-CoV-2 (2019-nCoV) Nucleocapsid
<b>Preservatives</b>	0.01M PBS(pH7.4) with 0.1% Proclin300
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	≥1 mg/ml
<b>Storage</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>MW</b>	46 kDa
<b>Uniprot ID</b>	<b>P0DTC2</b>
<b>Dilution Range</b>	WB=1:500-2000,ELISA=1:5000-10000
<b>Expiration Date</b>	12 months from date of receipt.