

Product Datasheet

Tick-borne encephalitis virus Genome polyprotein Antibody (orb1152260)

Description	Tick-borne encephalitis virus Genome polyprotein Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q01299
Expiration Date	12 months from date of receipt.