

Product Datasheet

Fission yeast BRL2 Antibody (orb1152295)

Description	Fission yeast BRL2 Antibody
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	O74563
Expiration Date	12 months from date of receipt.