



Product Datasheet

TRPV6 Antibody (orb1152309)

Description	TRPV6 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	E.coli-derived human TRPV6 recombinant protein (Position: A15-I765).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	83 kDa
Uniprot ID	Q9H1D0
Expiration Date	12 months from date of receipt.



Western blot analysis of TRPV6 using ant...



Flow Cytometry analysis of RT4 cells usi...