
Product Datasheet

**MPPB/PMPCB Antibody (monoclonal, 9F13E4)
(orb1152403)**

Description

MPPB/PMPCB Antibody (monoclonal, 9F13E4)


Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation


Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

E.coli-derived human MPPB/PMPCB recombinant protein (Position: E23-Q479).



IHC analysis of PMPCB using anti-PMPCB a...

Form/Appearance


Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.



IHC analysis of PMPCB using anti-PMPCB a...

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype


Mouse IgG2b

Clonality

Monoclonal

MW

43 kDa

Uniprot ID
075439


IF analysis of PMPCB using anti-PMPCB an...

Expiration Date

12 months from date of receipt.