

Product Datasheet

Drebrin/DBN1 Antibody (monoclonal, 4F6E7) (orb1152404)

Description

Drebrin/DBN1 Antibody (monoclonal, 4F6E7)

Species/Host

Mouse

Reactivity

Human, Monkey, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IHC, WB

Applications
Immunogen

E.coli-derived human Drebrin/DBN1 recombinant protein (Position: H9-D649).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2a

Clonality

Monoclonal

MW


120 kDa

Uniprot ID
Q16643
Sensitivity


> 5000 cells

Expiration Date


12 months from date of receipt.



IHC
analysis of
DBN1 using
anti-DBN1
ant...



IHC
analysis of
DBN1 using
anti-DBN1
ant...



IHC
analysis of
DBN1 using
anti-DBN1
ant...