
Product Datasheet

**Integrin alpha V/ITGAV Antibody (monoclonal,
8B10H2) (orb1152406)**

Description

Integrin alpha V/ITGAV Antibody (monoclonal,


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

ICC, IF, IHC, WB



IHC analysis of ITGAV using anti-ITGAV a...

Immunogen

E.coli-derived human Integrin alpha V/ITGAV recombinant protein (Position: H732-D970).

Form/Appearance


Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.



IHC analysis of ITGAV using anti-ITGAV a...

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2b

Clonality


Monoclonal

MW

130 kDa

Uniprot ID
P06756
Expiration Date

12 months from date of receipt.



IF analysis of ITGAV using anti-ITGAV an...