

Product Datasheet

NDFIP2 Antibody (monoclonal, 10D6D7) (orb1152407)

Description

NDFIP2 Antibody (monoclonal, 10D6D7)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

E.coli-derived human NDFIP2 recombinant protein (Position: M16-L336).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW


39 kDa

Uniprot ID


[Q9NV92](#)

Expiration Date


12 months from date of receipt.



IHC analysis of NDFIP2 using anti-NDFIP2...



IHC analysis of NDFIP2 using anti-NDFIP2...



IHC analysis of NDFIP2 using anti-NDFIP2...