


Product Datasheet

Human TFF2 ELISA Kit (DIY Antibody Pairs) (orb1152418)

Description	Human TFF2 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E24-Y129
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q03403
Expiration Date	Please enquire.



Human
TFF2 EZ
Set ELISA
Kit
standard
cur...