


Product Datasheet

Mouse IL-17 ELISA Kit (DIY Antibody Pairs) (orb1152426)

Description	Mouse IL-17 ELISA Kit (DIY Antibody Pairs)
Reactivity	Mouse
Range	15.6 pg/ml-1000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli, T22-A158
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q62386
Expiration Date	Please enquire.



Mouse IL-17 EZ-Set ELISA Kit standard CU...