

## Product Datasheet

### Rat B7-H4 Protein (orb1152428)

<b>Description</b>	Rat B7-H4, His Tag is expressed from human 293 cells (HEK293). It contains AA Phe 29 - Gly 257 (Acce...
<b>Reactivity</b>	Rat
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	B7-H4
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 27.1 kDa. The protein migrates as 45-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NM_001024244.1
<b>Purity</b>	90%
<b>Source</b>	Rat B7-H4, His Tag (orb1152428) is expressed from human 293 cells (HEK293). It contains AA Phe 29 - Gly 257 (Accession # Q501W4-1).
<b>MW</b>	27.1 kDa
<b>Uniprot ID</b>	<a href="#">Q501W4</a>
<b>NCBI</b>	<a href="#">NM_001024244.1</a>
<b>Expiration Date</b>	6 months from date of receipt.