


Product Datasheet

Cynomolgus IL-4 Protein (orb1152433)

Description	Cynomolgus IL-4, Fc Tag is expressed from human 293 cells (HEK293). It contains AA His 25 - Ser 153 ...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-4
Preservatives	51 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 41.4 kDa. The protein migrates as 46-50 kDa and 53 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Cynomolgus IL-4, Fc Tag (orb1152433) is expressed from human 293 cells (HEK293). It contains AA His 25 - Ser 153 (Accession # P79339-1).
MW	41.4 kDa
Uniprot ID	P79339
Expiration Date	6 months from date of receipt.



Immobilized
Cynomolgus
/ Rhesus
macaque ...