

Product Datasheet

Cynomolgus FGF R2 (IIIb) protein (orb1152442)

Description	Cynomolgus FGFR2 (IIIb), His Tag is expressed from human 293 cells (HEK293). It contains AA Pro 154 ...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	FGF R2 (IIIb)
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 25.9 kDa. The protein migrates as 42-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Pro 154 - Lys 368
Purity	95%
Source	Cynomolgus FGFR2 (IIIb), His Tag (orb1152442) is expressed from human 293 cells (HEK293). It contains AA Pro 154 - Lys 368 (Accession # A0A2K5TL84-1).
MW	25.9 kDa
Uniprot ID	A0A2K5TL84
Expiration Date	6 months from date of receipt.