

Product Datasheet

Mouse PVRIG Protein (orb1152453)

Description	Mouse PVRIG, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 35 - Asp 165 (Ac...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PVRIG
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.2 kDa. The protein migrates as 20-22 kDa and 23-28 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_011239268.1
Purity	90%
Source	Mouse PVRIG, His Tag (orb1152453) is expressed from human 293 cells (HEK293). It contains AA Ser 35 - Asp 165 (Accession # A0A1B0GS01-1).
MW	16.2 kDa
Uniprot ID	A0A1B0GS01
NCBI	XP_011239268.1
Expiration Date	6 months from date of receipt.