
Product Datasheet

RSV antibody (orb1152503)

Description

Mouse monoclonal antibody to RSV

Species/Host

Human

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA

Applications**Immunogen**

This antibody was prepared by the humanization of an F protein-specific murine MAb, RSV19. More specifically, molecular techniques were used to insert the complementarity determining regions from murine RSV19 into a human IgG1 heavy and k light chain-variable domain framework.

Target

RSV

Preservatives

PBS with 0.02% Proclin 300.

Concentration

1 mg/ml

Storage

Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.

Note

For research use only

Application notes

Though shown to be safe and well-tolerated (Everitt, 1996) and protective against RSV infection in BALB/c mice and cotton rats (Wyde, 1995), this antibody did not demonstrate efficacy in a Phase 3 study (Reichert, 2008). In in vitro studies, palivizumab, another humanized Mab specific for the F protein of RSV, was 4 to 5-fold more potent than felvizumab in neutralising RSV (Scott, 1999).

Isotype

IgG2a

Clonality

Monoclonal

Purity

Purified

Expiration Date

12 months from date of receipt.