

Product Datasheet

Integrin alpha 4 antibody (orb1152508)

Description	Mouse monoclonal antibody to Integrin alpha 4
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	The RAMOS cell line, a non-EBV releasing cell line, injected into mice.
Target	Integrin alpha 4
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds to human integrin alpha 4 and has been used to treat Multiple Sclerosis and Crohn's Disease.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P13612
Expiration Date	12 months from date of receipt.