

Product Datasheet

VEGF antibody (orb1152509)

Description	Mouse monoclonal antibody to VEGF
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC, IP, WB
Immunogen	Recombinant human VEGF.
Target	VEGF
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This is a research-grade biosimilar against VEGF. VEGF is a growth factor active in angiogenesis, vasculogenesis and endothelial cell growth. It induces endothelial cell proliferation, promotes cell migration, inhibits apoptosis and induces permeabilization of blood vessels.
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P15692
Expiration Date	12 months from date of receipt.