

Product Datasheet

TNF alpha antibody (orb1152525)

Description	Mouse monoclonal antibody to TNF alpha
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant human TNF (rTNF).
Target	TNF alpha
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody can be used to detect and neutralize human TNF alpha. Unlike Adalimumab (Ab00718), Infliximab does not bind to mouse TNF alpha (Arnoult 2017, PMID: 28584008).
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P01375
Expiration Date	12 months from date of receipt.