

## Product Datasheet

# Mouse Monoclonal anti-Human MICA/MICB Antibody (orb1154273)

<b>Description</b>	Mouse Monoclonal anti-Human MICA/MICB Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IHC, IP
<b>Immunogen</b>	Recombinant human MICA/MICB
<b>Storage</b>	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG2a
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.