

Product Datasheet

Mouse Monoclonal anti-Human CD273 (B7-DC, PD-L2) Antibody (orb1154516)

Description	Mouse Monoclonal anti-Human CD273 (B7-DC, PD-L2) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	Human PD-L2 cDNA
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.