

Product Datasheet

Mouse Monoclonal anti-His Tag Antibody (orb1154688)

Description	Mouse Monoclonal anti-His Tag Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	KLH-conjugated His-tag peptide.
Storage	Reconstitute the antibody with 200 ul sterile PBS and the final concentration is 0.5mg/ml. Reconstituted antibody can also be aliquoted and stored frozen at -2°C for at least for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.