

## Product Datasheet

# Goat anti-Mouse IgG(H+L), Cross absorbed against human and rat immunoglobulins, Purified (orb1154766)

<b>Description</b>	Goat anti-Mouse IgG(H+L), Cross absorbed against human and rat immunoglobulins, Purified
<b>Conjugation</b>	Unconjugated
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.