

Product Datasheet

Goat anti-Mouse IgG(H+L), Cross absorbed against human and rabbit immunoglobulins, Purified (orb1154767)

Description	Goat anti-Mouse IgG(H+L), Cross absorbed against human and rabbit immunoglobulins, Purified
Conjugation	Unconjugated
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.