

Product Datasheet

Rabbit anti FLAG tag (DYKDDDDK) Polyclonal Antibody (orb1154822)

Description	Rabbit anti FLAG tag (DYKDDDDK) Polyclonal Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	FLAG-tag
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.