

Product Datasheet

Mouse anti HSV-tag Monoclonal Antibody (orb1154828)

Description	Mouse anti HSV-tag Monoclonal Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	HSV-tag
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.