

Product Datasheet

Glu-Glu Tag Mouse Monoclonal Antibody(2G1) (orb1154835)

Description	Glu-Glu Tag Mouse Monoclonal Antibody(2G1)
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.