

Product Datasheet

Rsk-1/2/3/4 (phospho Ser221/227/S218/232) Polyclonal Antibody (orb1155817)

Description	Rsk-1/2/3/4 (phospho Ser221/227/S218/232) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.