

Product Datasheet

MLK1/2 (phospho Thr312/266) Polyclonal Antibody (orb1156455)

| | |
|----------------------------|---|
| Description | MLK1/2 (phospho Thr312/266) Polyclonal Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC-P |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Expiration Date | 12 months from date of receipt. |