

Product Datasheet

Flk-1/Flt-4 (phospho Tyr1054/Y1063) Polyclonal Antibody (orb1156629)

Description	Flk-1/Flt-4 (phospho Tyr1054/Y1063) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.