

Product Datasheet

SLC27A2 antibody (orb11575)

Description Rabbit polyclonal antibody to SLC27A2

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen Recombinant fusion protein of human FATP2(NP_003636.2).

Preservatives PBS with 0.05% proclin300,50% glycerol,pH7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 70KDa


Uniprot ID [O14975](#)

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of human liver tissue.



Western blot analysis of extracts of HepG2 cells.