

Product Datasheet

Na⁺/K⁺-ATPase α 1 (Phospho-Tyr260) Polyclonal Antibody (orb1157656)

| | |
|----------------------------|---|
| Description | Na ⁺ /K ⁺ -ATPase α 1 (Phospho-Tyr260) Polyclonal Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Expiration Date | 12 months from date of receipt. |