

Product Datasheet

MER/TYRO3 (Phospho-Tyr753/Tyr685) Polyclonal Antibody (orb1157855)

Description	MER/TYRO3 (Phospho-Tyr753/Tyr685) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.