

Product Datasheet

Phospho-ERK 1/2 (Y222/205) Mouse Monoclonal Antibody(4G3) (orb1161149)

Description	Phospho-ERK 1/2 (Y222/205) Mouse Monoclonal Antibody(4G3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.