

Product Datasheet

Phospho-Akt (S473) Mouse Monoclonal Antibody(6F8) (orb1162360)

Description	Phospho-Akt (S473) Mouse Monoclonal Antibody(6F8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.