

Product Datasheet

V5-Tag Monoclonal Antibody (orb1162600)

| | |
|----------------------------|---|
| Description | V5-Tag Monoclonal Antibody |
| Species/Host | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IF, IP, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |