

Product Datasheet

Splicing factor 1 (Phospho-Ser82) Polyclonal Antibody (orb1164510)

Description	Splicing factor 1 (Phospho-Ser82) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.