

Product Datasheet

p70 S6 kinase α (Phospho-Ser418) Polyclonal Antibody (orb1164512)

Description	p70 S6 kinase α (Phospho-Ser418) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.