

Product Datasheet

ATP1 alpha1/Na+K+ ATPase1 (Phospho-Ser23) Antibody (orb1164775)

Description	ATP1 alpha1/Na+K+ ATPase1 (Phospho-Ser23) Antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.