

Product Datasheet

CDK1/CDC2 (Phospho-Thr14) Antibody (orb1164830)

Description CDK1/CDC2 (Phospho-Thr14) Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Clonality Polyclonal

Expiration Date 12 months from date of receipt.