

Product Datasheet

Ephrin B1/B2/B3 (Phospho-Tyr324) Antibody (orb1164985)

Description	Ephrin B1/B2/B3 (Phospho-Tyr324) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.