

Product Datasheet

RSK1/2/3/4 (Phospho-Ser221/227/S218/232) Antibody (orb1165011)

Description	RSK1/2/3/4 (Phospho-Ser221/227/S218/232) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.