

Product Datasheet

Interferon- α Receptor Type I Subunit I (phospho-Ser535/Ser539) Antibody (orb1165129)

Description	Interferon- α Receptor Type I Subunit I (phospho-Ser535/Ser539) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.