

Product Datasheet

AMPK α 1/AMPK α 2(Phospho-Thr183/Thr172) Antibody (orb1165291)

Description	AMPK α 1/AMPK α 2(Phospho-Thr183/Thr172) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.