

## Product Datasheet

# Estrogen Receptor- $\alpha$ (Phospho-Ser167) Antibody (orb1165393)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Estrogen Receptor- $\alpha$ (Phospho-Ser167) Antibody                                                                                   |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Reactivity</b>          | Human                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | IF, IHC, WB                                                                                                                             |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |